

Supporting Information

Lin et al. 10.1073/pnas.1115025109




Fig. S1. SopA/Ubch7 pulldown experiment. His-SopA was incubated with Ubch7 in a 1:1 molar ratio for 30 min at 4 °C. Co^{2+} beads were then added to the mixture and incubated for 10 min before centrifugation. After extensive washing, proteins bound to the beads were analyzed on SDS-PAGE with Coomassie blue staining and shown in Fig. 2B in the main text. The input samples and unbound fractions are shown here per request of a reviewer.

