



Supplemental Figure S4. Effect of DHC on early and late stage of HUVECs apoptosis was detected by flow cytometry with Annexin-V-FITC/PI dual staining. A representative histogram of flow cytometric analysis using double staining with annexin-V (FITC-A) and PI (PE-A). HUVECs were treated with DHC (5 μM) in EBM-2 basal medium for 4hr or 24hr. The lower right quadrants represent the cells in the early stage of apoptosis. The upper right plus left quadrants contain the cells in the late stage of apoptosis and necrosis. Data represent from three independent experiments.