

Table S1: RT-PCR primers used for MuLV analysis		
Gene	Forward primer	Reverse primer
MuLV	5'-CTGAGAGAAGTCAACAAGCG-3'	5'-CTGGATCTCTCCACTCAAAGGC-3'
GAPDH	5'-GAAATCCCATCACCATCTTCCAGG-3'	5'-GAGCCCCAGCCTTCTCCATG-3'
<u>MuLV Junction</u>	Forward Primer	Reverse Primer
Junc 6	5'-GCTCGTCCGGGATTTGGAGA-3'	5'-CAGCAACTGTTGACAGTCATCCCA-3'
Junc 8	5'-GCTCGTCCGGGATTTGGAGA-3'	5'-TTCCAGACCTCTCGTTGTACTGCT-3'
Junc 128	5'-AATTGCCCAATGAGATCGAGGCTG-3'	5'-TTCCAGACCTCTCGTTGTACTGCT-3'
Junc 225	5'-ATGACTTACTGCTGGCCGCTACTT-3'	5'-GGTTCCATGCGTTTCAGCCAAGAT-3'
Reference for GAPDH primer: West et al, JBC 2004, Vol 279, (28): 28896-28902		