

Supplemental Table

Primers used in real-time quantitative PCR

Human 2',5'-Oligoadenylate synthase-1A	Sense	5' -GTGTGTCCAAGGTGGTAAAGG-3'
	Antisense	5' -CTGCTCAAACCTTCACGGAA-3'
Human Mx1	Sense	5' -GCCAGGACCAGGTATACAG-3'
	Antisense	5' -GCTCCTTCAGGAGCCAGA -3'
Human IRS-2	Sense	5' -AGCCGGGAGACAAGAG-3'
	Antisense	5' -TGATACTCTCTGACATGT-3'
Human GAPDH	Sense	5' -GCCAAAAGGGTCATCATCT-3'
	Antisense	5' -ATGGATGACCTTGCCAG-3'
Mouse 2',5'-Oligoadenylate synthase-1A	Sense	5' -GCTGATGTCAAATCAGCCGT-3'
	Antisense	5' -AGCTTGAAGCTCAGAGACC-3'
Mouse Mx1	Sense	5' -GGGGAGGAAATAGAGAAAATGAT-3'
	Antisense	5' -GTTTACAAAGGGCTTGCTTGCT-3'
Mouse PTP1B	Sense	5' -CGGCTATTTACCAGGACATTC-3'
	Antisense	5' -TGCGGTTGAGCATGACCAC-3'
Mouse β-Actin	Sense	5' -ATGGATGACGATATCGCTG-3'
	Antisense	5' -GTTGGTAACAATGCCATGTTC-3'