



Figure S4 Physical interaction between Endos and B55/Twins *in vitro*. Pull-down experiment using glutathione-Sepharose beads bound with the same GST fusion proteins as in Figure 6, but incubated with recombinant, purified Twins protein instead of *in vitro* transcribed/translated ³⁵S-labeled Twins. (Top) Twins protein binding to beads, as detected on Western blots (WB). Depending on the sample, between 10% and 20% of the input Twins was retained on the beads (not shown). (Bottom) Coomassie blue-stained gel showing that the beads used for *in vitro* binding reactions had similar amounts of bound GST::Endos fusion proteins.