

Table B. Primers used in quantitative polymerase chain reactions.

Gene Product	Forward Primer	Reverse Primer
GSK3 β	5'-AGAGATAAAGATGGCAGCAAGG-3'	5'-ATCCATTTCCAATCACTTTTCG-3'
GSK3 α	5'-ACAGCGGGAAGGTGACCACAGTG-3'	5' AAGTACCGCAGCCTCACAATA-3'
GSTM5	5'-GAGAAGCAGAAGCCAGAGT-3'	5'-ATACGATACTGGTCAAGAAT-3'
HO-1	5'-CACAGATGGCGTCACTTCGTC-3'	5'-GTGAGGACCCACTGGAGGAG-3'
GCLM	5'-TTACCGAGGCTACGTGTCAGAC-3'	5'-TATCGATGGTCAGGTCGATGTC-3'
GCLC	5'-AATCAGCCCCGATTTAGTCAGG-3'	5'-CCAGCTGTGCAACTCCAAGGAC-3'
β -actin	5'-TCCTTCCTGGGCATGGAG-3'	5'-AGGAGGAGCAATGATCTTGATCTT-3'