

Table S1

Supplemental Table S1. Primers used in conventional PCR experiments

GENE	Sense Primer	Anti-Sense Primer	Expected Size (bps)
UGT2B4	5'-CTCTCCTGGCTACGCAA-3'	5'-TTTGCAGTGGAGTCCTCC-3'	326
UGT2B7	5'-CTTCTCTCCTGGCTACACTT-3'	5'-TTTGCAGTGGAGTCCTCC-3'	329
UGT2B10	5'-CTGGCTGAGCTATTTAACATACCC-3'	5'-CGCATAAGCCATATGTCAGC-3'	298
UGT2B11	5'-GACATCGTTTTTGCAGATGCTGT-3'	5'-CGCATAAGCCATATGTCAGC-3'	306
UGT2B15	5'-GTTGGGAATATTATGACTACAGTAAC-3'	5'-TTTGGTAAGAATGGGCGAG-3'	462
UGT2B17	5'-GTGTTGGGAATATTCTGACTATAATATA-3'	5'-TTTGGTAAGAATGGGCGAG-3'	462
UGT2B28	5'-CAGTCTCTGCTTCACTCC-3'	5'-TTTGCAGTGGAGTCCTCC-3'	322
UGT1A	5'-TCCTGTGGCGGTACACTG-3'	5'-TCAATGGGTCTTGGATTTGTGG-3'	600
GAPDH	5'-CCAGGGCTGCTTTTAACTC-3'	5'-GCTCCCCCTGCAAATGA-3'	291