

Supplementary information, Table S4 Primer sequences (5' → 3') for ChIP assay

Name	5' primers	3' primers
<i>Id1</i>	TTCAGACGCTGACACAGACC	CCTCTTGCAGAAAGGTTCCA
<i>Msx2</i>	CCGAGTATCTACCTAAATTCCTGCTG	TTCCTTCTAATGGGCCGCTTGTT
<i>Id1</i> -3'UTR	TTGGGAAAGCGGGTCACTGGA	TGGCTCAGCCGGAAGAGGG
<i>AP2γ</i> -P1	GGGTAGGGCAGTCTCAGGATT	GAACCCTTGACTTTACCAGGCTC
<i>AP2γ</i> -P2	CTGCGCTTAGTTCTGGAAGCTC	AATGGGCAAGCACCGGACTTTG
<i>AP2γ</i> -P3	CAGGCTGCTAGGTAGTTGACTG	CACCCAAGGCCAAAGAAGTCAGT
<i>AP2γ</i> -P4	TGTGGGTACCCCAAGCCTAAG	CACTGCTTCTAACTCACCCAC