

SUPPLEMENTAL MATERIAL

Zheng et al., <http://www.jem.org/cgi/content/full/jem.20120342/DC1>**Table S1.** ChIP Assay primers used to confirm Egr2 association

Gene	Primers
<i>DGK-ζ</i>	Forward, 5'-GCAGACCTAGGACACGGTAGA-3'; Reverse, 5'-CTGTGGGCAGAAACAAACCT-3'
<i>Cbl-b</i>	Forward, 5'-TAGGAAGGTGAGTGAGGCAGGAAG-3'; Reverse, 5'-CCAATGGGCGGGACGAGAC-3'
<i>Itch</i>	Forward, 5'-AGTGTGGCGCTCGAGTCT-3'; Reverse, 5'-GCCAAGCACTAGCAGTTCG-3'
<i>Tob1</i>	Forward, 5'-GGGAGGAGGCTATAAAACGAG-3'; Reverse, 5'-GCGCTCGCTTCTTACACA-3'
<i>Dtx1</i>	Forward, 5'-CCTCCTGACACCCATCTC-3'; Reverse, 5'-GGGATTTCCTGGAGTTCAT-3'

Table S2. qRT-PCR primers and probes used to confirm dependency on Egr2 for expression in anergy

Gene	Primers	Probe
<i>DGK-ζ</i>	Forward, 5'-CCCAAGGTGAAGAACTGATTG-3'; Reverse, 5'-TCTTCTGGTCACGGTGCAT-3'	5'-GAGCAGGA-3' (Roche)
<i>Cbl-b</i>	Forward, 5'-ATGTCCTCCTCGGCTTT-3'; Reverse, 5'-GAGACAATTTGCTAATGGACCAG-3'	5'-CTTCTCCC-3' (Roche)
<i>Itch</i>	Ordered from Applied Biosystems	Ordered from Applied Biosystems
<i>Tob1</i>	Forward, 5'-TCCAGCTCTCCCTCACCTC-3'; Reverse, 5'-ACCGGGCATGAAGGTAG-3'	5'-TCTGCTGC-3' (Roche)
<i>Dtx1</i>	Forward, 5'-CCAACCCAGGCAAGAGTT-3'; Reverse, 5'-GATGAGCAATCTCAGCACCTT-3'	5'-AGCCAGAG-3' (Roche)