


A new, rapid and reproducible method to obtain high quality endothelium *in vitro*

Cytotechnology

Nuria Jiménez, Vincent J.D. Krouwer and Jan A. Post

HUVECs-p1: 7-days cobblestone



Online Resource 2 *Phase contrast microscopy of 7-days cobblestone HUVECs on Matrigel, fibronectin and gelatin.* Cells (passage 1, p1) were seeded at 20000 cells/cm² on Aclar coated with different matrices and left to reach the 7-days cobblestone state. In all the cases, cells show distinct dark cell limits associated to the mature cobblestone appearance. Scale bar (applicable to all the panels): 50 μ m