

Additional file 2 Table of primers used for RT-PCR analysis

Primer name	Sequence	Amplicon length
<i>fkB-1</i> _5	GCGGCGGCGTTGGCACCTCAT	280
<i>fkB-1</i> _3	ACCCGCGCCCCGTCTCCAAACTC	
<i>fkBG</i> _5	GTACGGATCGGCGACGCCCTGGAC	258
<i>fkBG</i> _3	GTCGCGCAGCCCCGCATTGAGTG	
<i>fkBN</i> _5	GGTGAACTCGGCCCGCAGC	339
<i>fkBN</i> _3	CGAGCAGCCCCTTATGAACTGATCGC	
<i>allA-1</i> _5	GCCCTGGCGGAGCTTCTGGAGTCA	283
<i>allA-1</i> _3	TCACCGTCTGGCTGGGGGAGTTCA	
<i>allA-2</i> _5	GGCGGGATGCTGGCGGTGAAGG	365
<i>allA-2</i> _3	CGGGTCATGTCGTTCCACCATCAGG	
<i>fkBR</i> _5	GCCGACATCGCCCTCCAGGTG	257
<i>fkBR</i> _3	GGGTCCTGCGGGGCCTCTGC	
<i>hrdB</i> _5	GCCGCCGCGCCAAGAACC	419
<i>hrdB</i> _3	CCAGCGGCGTGTGCAGCGAGA	