

**Table S3**

TaqMan® MicroRNA Assays Human Panel - Early Access Kit		Species Targeted by Assay (Sanger miRNA Registry 9/04, Version 5.0)		
AB Assay Name	Mature miRNA Sequence (Assay Target)	Human	Mouse	Rat
hsa-let-7a	ugagguaguagguguauaguuu	Yes	Yes	Yes
hsa-let-7b	ugagguaguaggugugugguu	Yes	Yes	Yes
hsa-let-7d	agagguaguaggugcauagu	Yes	Yes	Yes
hsa-let-7e	ugagguaggagguguauagu	Yes	Yes	Yes
hsa-let-7g	ugagguaguaguuuguacagu	Yes	Yes	
hsa-let-7i	ugagguaguaguuugugcu	Yes	Yes	Yes
hsa-miR-9	ucuuuggguuaucuagcguuuga	Yes	Yes	Yes
hsa-miR-9*	uaaaagcuaguaaacccgaaagu	Yes	Yes	
hsa-miR-10a	uaccgcguagauccggaaauuugug	Yes	Yes	Yes
hsa-miR-15a	uagcagcacauuaugguuugug	Yes	Yes	
hsa-miR-15b	uagcagcacaucaugguuuaca	Yes	Yes	Yes
hsa-miR-16	uagcagcacguaaaaauuuggcg	Yes	Yes	Yes
hsa-miR-17-3p	acugcagugaaggcacuugu	Yes		
hsa-miR-17-5p	caaagugcuuacagugcaggagu	Yes	Yes	Yes
hsa-miR-19a	ugugcaaaaucuaugcaaaacuga	Yes	Yes	Yes
hsa-miR-20	uaaagugcuuuaugugcaggua	Yes		
hsa-miR-21	uagcuuauacagacugauguuga	Yes	Yes	Yes
hsa-miR-23a	aucacauuugccagggauuucc	Yes	Yes	Yes
hsa-miR-23b	aucacauuugccagggauuaccac	Yes	Yes	Yes
hsa-miR-25	cauugcacuugucucggucuga	Yes	Yes	Yes
hsa-miR-26a	uucaaguuaauccaggauaggcu	Yes	Yes	Yes
hsa-miR-26b	uucaaguuaauucaggauagg	Yes		
hsa-miR-27a	uucacaguggcuaaguuccgc	Yes		
hsa-miR-27b	uucacaguggcuaaguucug	Yes	Yes	Yes
hsa-miR-28	aaggagcucacagucuaauugag	Yes	Yes	Yes
hsa-miR-29a	cuagcaccaucugaaaucgguu	Yes	Yes	Yes
hsa-miR-29b	uagcaccuuuugaaaucagu	Yes		
hsa-miR-29c	uagcaccuuuugaaaucggua	Yes	Yes	Yes
hsa-miR-30a-3p	cuuucagucggauuuugcagc	Yes	Yes	Yes
hsa-miR-30b	uguaaacauccuacacucagc	Yes	Yes	Yes
hsa-miR-30c	uguaaacauccuacacucucagc	Yes	Yes	Yes
hsa-miR-30d	uguaaacauccccgacuggaag	Yes	Yes	Yes
hsa-miR-30e	uguaaacauccuugacugga	Yes	Yes	Yes
hsa-miR-31	ggcaagaugcuggcauagcug	Yes		
hsa-miR-34a	uggcagugucuuuagcugguugu	Yes		
hsa-miR-34b	aggcagugucuuuagcugauug	Yes		
hsa-miR-34c	aggcaguguguuuagcugauug	Yes		
hsa-miR-92	uuuugcacuuuguccggccugu	Yes		
hsa-miR-95	uucaacggguauuuauugagca	Yes		
hsa-miR-96	uuuggcacuagcacaauuuuugc	Yes		
hsa-miR-98	ugagguaguaaguuguauuuuu	Yes	Yes	Yes
hsa-miR-99a	aacccguagauccgaucuugug	Yes		Yes
hsa-miR-100	aacccguagauccgaacuugug	Yes	Yes	Yes
hsa-miR-103	agcagcauuguacagggcuauga	Yes	Yes	Yes
hsa-miR-104	ucaacaucagucugauaagcua	Yes		

hsa-miR-105	ucaaaugcucagacuccugu	Yes		
hsa-miR-106a	aaaagugcuuacagugcaggguagc	Yes		
hsa-miR-107	agcagcauuguacagggcuaauca	Yes	Yes	Yes
hsa-miR-122a	uggagugugacaauuggguuuugu	Yes	Yes	Yes
hsa-miR-124a	uuuaggcacgcggugaaugccaa	Yes	Yes	Yes
hsa-miR-124b	uuuaggcacgcggugaaugc	Yes		
hsa-miR-125a	ucccugagaccuuuaaccugug	Yes	Yes	Yes
hsa-miR-125b	ucccugagaccuuacuuuguga	Yes	Yes	Yes
hsa-miR-126	ucguaccgugaguauuaaugc	Yes	Yes	Yes
hsa-miR-127	ucggaucgcugagcugcuuggcu	Yes	Yes	Yes
hsa-miR-128a	ucacagugaaccggucucuuuu	Yes	Yes	Yes
hsa-miR-128b	ucacagugaaccggucucuuuc	Yes	Yes	Yes
hsa-miR-129	cuuuuugcggucugggcuugc	Yes		
hsa-miR-130a	cagugcaauguuuaagggc	Yes	Yes	Yes
hsa-miR-130b	cagugcaaugaugaaagggcau	Yes	Yes	Yes
hsa-miR-132	uaacagucuacagcccaugggucg	Yes	Yes	Yes
hsa-miR-133a	uuggccccuuaaccaggcugu	Yes	Yes	Yes
hsa-miR-133b	uuggccccuuaaccaggcua	Yes	Yes	
hsa-miR-134	ugugacugguugaccagaggg	Yes	Yes	Yes
hsa-miR-135a	uauggcuuuuuauuccuauguga	Yes	Yes	Yes
hsa-miR-135b	uauggcuuuuucauuccuaugug	Yes	Yes	Yes
hsa-miR-137	uaauugcuaaagaauacgcguag	Yes	Yes	Yes
hsa-miR-138	agcugguguguguaauc	Yes	Yes	Yes
hsa-miR-139	ucuacagugcacgugucu	Yes	Yes	Yes
hsa-miR-140	agugguuuuuacccuaugguag	Yes	Yes	Yes
hsa-miR-141	aacacugucugguaaaagaugg	Yes	Yes	
hsa-miR-142-3p	uguaguguuuccuacuuuaugga	Yes		
hsa-miR-142-5p	cauaaaguagaaagcacuac	Yes	Yes	Yes
hsa-miR-144	uacaguauagaugauaguacuag	Yes	Yes	Yes
hsa-miR-145	guccaguuuuucccaggaauccuu	Yes	Yes	Yes
hsa-miR-146	ugagaacugaaauuccaugggu	Yes	Yes	Yes
hsa-miR-147	guguguggaaaugcuucugc	Yes		
hsa-miR-148a	ucagugcacuacagaacuuugu	Yes	Yes	
hsa-miR-149	ucuggcuccgugucuucacucc	Yes	Yes	
hsa-miR-150	ucucccaacccuuguaccagug	Yes	Yes	Yes
hsa-miR-151	acuagacugaaggcuccuugagg	Yes		
hsa-miR-152	ucagugcaugacagaacuugg	Yes	Yes	Yes
hsa-miR-154	uagguaauccguguugccuucg	Yes	Yes	Yes
hsa-miR-154*	aaucauacacgguugaccuauu	Yes		
hsa-miR-155	uuuaugcuaaucugugauagggg	Yes		
hsa-miR-181a	aacauucaacgcugucggugagu	Yes	Yes	Yes
hsa-miR-181b	aacauucauugcugucgguggguu	Yes	Yes	Yes
hsa-miR-181c	aacauucaaccugucggugagu	Yes	Yes	Yes
hsa-miR-182	uuuggccaugguagaacucaca	Yes	Yes	
hsa-miR-182*	ugguucuagacuugccaaacua	Yes		
hsa-miR-183	uauggcacugguagaauucacug	Yes	Yes	Yes
hsa-miR-184	uggacggagaacugauaagggu	Yes	Yes	Yes
hsa-miR-185	uggagagaaaggcaguuc	Yes	Yes	Yes
hsa-miR-186	caaagaauucuccuuuuggcuu	Yes	Yes	Yes
hsa-miR-187	ucgugucuuguguugcagccg	Yes		
hsa-miR-189	gugccuacugagcugauaucagu	Yes	Yes	

hsa-miR-190	ugauauguuugauauuuagggu	Yes	Yes	Yes
hsa-miR-191	caacggaauccaaaggcagcu	Yes	Yes	Yes
hsa-miR-193	aacuggccuacaaaagucccag	Yes	Yes	Yes
hsa-miR-194	uguaacagcaacuccauggga	Yes	Yes	Yes
hsa-miR-195	uagcagcacagaaaauuuggc	Yes	Yes	Yes
hsa-miR-197	uucaccaccuucuccaccaggc	Yes		
hsa-miR-198	gguccagaggggagauagg	Yes		
hsa-miR-199a	cccaguguucagacuaccuguuc	Yes	Yes	Yes
hsa-miR-199a*	uacaguagucugcacauugguu	Yes	Yes	
hsa-miR-199b	cccaguguuuagacuaucuguuc	Yes		
hsa-miR-199-s	cccaguguucagacuaccuguu	Yes		
hsa-miR-200a	uaacacugucugguaacgaugu	Yes	Yes	Yes
hsa-miR-200b	cucuaauacugccugguaaugaug	Yes		Yes
hsa-miR-200c	aaucacugccgguaaugaugga	Yes	Yes	Yes
hsa-miR-203	gugaaauguuuaggaccacuag	Yes		Yes
hsa-miR-204	uucccuuugucauccuaugccu	Yes		Yes
hsa-miR-205	uccuucauuccaccggagucug	Yes	Yes	Yes
hsa-miR-210	cugugcgugugacagcggcug	Yes	Yes	Yes
hsa-miR-211	uucccuuugucauccuucgcu	Yes		
hsa-miR-213	accaucgaccguugauuguacc	Yes	Yes	Yes
hsa-miR-214	acagcaggcacagacaggcag	Yes	Yes	Yes
hsa-miR-215	augaccuaugaauugacagac	Yes		
hsa-miR-216	uaaucucagcuggcaacugug	Yes	Yes	Yes
hsa-miR-218	uugugcuuugaucuaaccaugu	Yes	Yes	Yes
hsa-miR-219	ugauuguccaaacgc当地	Yes	Yes	Yes
hsa-miR-220	ccacaccguaucugacacuuu	Yes		
hsa-miR-221	agcuacauugucugcugguuuc	Yes		Yes
hsa-miR-222	agcuacaucuggcuacugggucuc	Yes	Yes	Yes
hsa-miR-223	ugucaguuugucaaauacccc	Yes	Yes	Yes
hsa-miR-224	caagucacuagugguuccguuua	Yes		
hsa-miR-296	agggcccccccucaauccugu	Yes	Yes	Yes
hsa-miR-299	ugguuuaccgucccacauacau	Yes	Yes	Yes
hsa-miR-301	cagugcaauaguauugucaaagc	Yes	Yes	
hsa-miR-302a	uaagugcuuccauguuuugguga	Yes	Yes	
hsa-miR-302b	uaagugcuuccauguuuugguag	Yes		
hsa-miR-302b*	acuuuaacauggaagugc当地	Yes		
hsa-miR-302c	uaagugcuuccauguuucagugg	Yes		
hsa-miR-302c*	uuuaacauggggguaccugcug	Yes		
hsa-miR-302d	uaagugcuuccauguuugagugu	Yes		
hsa-miR-320	aaaagcuggguugagagggcgaa	Yes	Yes	Yes
hsa-miR-323	gcacauuacacggcugaccucu	Yes	Yes	Yes
hsa-miR-324-5p	cgcauccccuagggcauuggugu	Yes	Yes	Yes
hsa-miR-325	ccuaguagguguccaguaagu	Yes		
hsa-miR-326	ccucuggcccuuccuccag	Yes		
hsa-miR-328	cuggccucucugccuuccgu	Yes	Yes	Yes
hsa-miR-330	gcaaaggcacacggccugcagaga	Yes		
hsa-miR-331	gccccuggccuauccuugaa	Yes	Yes	Yes
hsa-miR-335	ucaagagcaauaacgaaaaaugu	Yes	Yes	Yes
hsa-miR-337	uccagcuccuaauaugaugccuuu	Yes		
hsa-miR-338	uccagcaucagugauuuuguuga	Yes	Yes	Yes
hsa-miR-339	ucccuguccuccaggagcuca	Yes	Yes	Yes

hsa-miR-340	uccgucucaguuacuuuuauagcc	Yes	Yes	Yes
hsa-miR-342	ucucacacagaaaucgcacccguc	Yes	Yes	Yes
hsa-miR-367	aauugcacuuuagcaaugguga	Yes		
hsa-miR-368	acauagaggaaauuccacguuu	Yes		
hsa-miR-370	gccugcuggguggaaccugg	Yes		
hsa-miR-371	gugccccaucuuuugagugu	Yes		
hsa-miR-372	aaagugcugcgacauuugagcgu	Yes		
hsa-miR-373	gaagugcuucgauuuugggugu	Yes		
hsa-miR-373*	acucaaaaugggggcguuucc	Yes		
hsa-miR-374	uuauuaauacaaccugauaagug	Yes		
ath-miR159a	uuuggauugaagggagcucua	<i>Arabidopsis</i> negative control		
cel-lin-4	ucccugagaccucaaguguga	<i>C. elegans</i> negative control		
cel-miR-2	uaucacagccagcuuugaugugc	<i>C. elegans</i> negative control		