



AM2



AM1

Supplementary information, Figure S2 Yeast-two-hybrid assay indicates that AM1 interact with almost all ABA receptors.

Interactions of 14 PYR/PYLs fused with binding-domain (BD) with activation-domain (AD)-fused HAB1 on -Leu/-Trp/-His SD media containing (+)-ABA and corresponding compounds were determined by yeast-two-hybrid assay. Working concentration is 2 μ M for AM1, AM2 and (+)-ABA. -Leu/-Trp (-L/-T) and -Leu/-Trp/-His SD media containing 0.1%DMSO (DMSO) are used as positive and negative controls, respectively.