

Multi-channel PMMA microfluidic biosensor with integrated IDUAs for electrochemical detection

Analytical and Bioanalytical Chemistry

Electronic Supplementary Material

Nongnoot Wongkaew, Peng He, Vanessa Kurth, Werarak Surareungchai, Antje J. Baeumner*

* Corresponding author, Department of Biological and Environmental Engineering, Cornell University, 202 Riley Robb Hall, Ithaca, NY 14853, USA, e-mail: ajb23@cornell.edu

Table S1 Effect of bubble generation during surface modification on IDUA










Method	PMMA surface activation by coupling reagents	Cystamine conjugation	IDUAs after etching
1	 2 mL of 300 mM Sulfo-NHS + 300 mM EDC in 0.05 M MES buffer	 2 mL of 300 mM cystamine in 0.05 M Carbonate buffer	 100x short circuited
2	 1 mL of 300 mM Sulfo-NHS + 300 mM EDC + 0.1% (v/v) Tween 20 in 0.05 M MES buffer	 1 mL of 300 mM cystamine + 0.1% (v/v) Tween 20 in 0.05 M Carbonate buffer	 100x intact circuit
3	 2 mL of 300 mM NHS + 300 mM EDC in 0.05 M MES buffer	 2 mL of 300 mM cystamine in 0.05 M Carbonate buffer	 100x short circuited

Table S1 Effect of bubble generation during surface modification on IDUA (continued)







Method	PMMA surface activation by coupling reagents	Cystamine conjugation	IDUAs after etching
5	 1 mL of 300 mM NHS + 300 mM EDC + 0.1% (v/v) Tween 20 in 0.05 M MES buffer	 1 mL of 300 mM cystamine + 0.1% (v/v) Tween 20 in 0.05 M Carbonate buffer	 short circuited
6	 30 mL of 300 mM NHS + 300 mM EDC in 0.05 M MES buffer	 30 mL of 300 mM cystamine in 0.05 M Carbonate buffer	 short circuited

Table S2 Optimization of bonding condition in which 4000 lbf bonding force, 8 min UV treatment time, 5 min bonding time with photoresist protection were used in the end. The first column provides conditions used. In the upper images for each condition a large view of the IDUAs with overlaid PMMA channels are shown. In the lower images a close-up of the IDUAs are shown. Defects are visible for all of the non-optimal conditions.

Parameter	Result	
<p>Bonding pressure</p> <ul style="list-style-type: none"> - UV treatment for 8 min - bonding time for 5 min - Without photoresist protection 	<p style="text-align: center;">4000 lbf</p>	<p style="text-align: center;">5000 lbf</p>
<p>Effect of photoresist protection</p> <ul style="list-style-type: none"> - UV treatment for 8 min - bonding time for 5 min - Applied bonding pressure 4000 lbf 	<p style="text-align: center;">With photoresist</p>	<p style="text-align: center;">Without photoresist</p>
<p>UV treatment time</p> <ul style="list-style-type: none"> - bonding time for 5 min - Applied bonding force 4000 lbf - Without photoresist protection 	<p style="text-align: center;">4 min</p>	<p style="text-align: center;">8 min</p>