

Supporting Table S2. Experimental conditions

Experimental conditions; Zeiss LSM 510 Meta Live

37 °C, 5% CO₂

Objective	63X/1.4 OIL PLAN APO
Scan mode	xyz
Format	256 x 256
ROI	2.5 µm in diameter
Argon laser power	50 mW, 100% output
Bleaching 488 nm	100% (Zoom in mode)
Imaging 488 nm	5%
Emission range	490-530 nm
Pre-bleach	10 sec
Bleach	7.5 sec
Post-bleach 1	80sec