

Table S2. Setup of the DIGE experiments. Letters A–D refer to biological replicate samples. Gels 1–4, 5–8 and 9–12 each form a separate DIGE experiment.

Gel no.	Cy3	Cy5	Cy2
1	4h A	20h A	Mix (4h A–D + 20h A–D)
2	4h B	20h B	Mix (4h A–D + 20h A–D)
3	20h C	4h C	Mix (4h A–D + 20h A–D)
4	20h D	4h D	Mix (4h A–D + 20h A–D)
5	12h A	20h A	Mix (12h A–D + 20h A–D)
6	12h B	20h B	Mix (12h A–D + 20h A–D)
7	20h C	12h C	Mix (12h A–D + 20h A–D)
8	20h D	12h D	Mix (12h A–D + 20h A–D)
9	20h A	28h A	Mix (20h A–D + 28h A–D)
10	20h B	28h B	Mix (20h A–D + 28h A–D)
11	28h C	20h C	Mix (20h A–D + 28h A–D)
12	28h D	20h D	Mix (20h A–D + 28h A–D)