

Electronic Supplementary Material

Supplementary Results

1. Bland Altman plot for Roche v Mercodia assay is shown in supplementary figure 1.

2. Detection of low concentrations of C-peptide by the Roche and Mercodia assay.

In the 67 samples selected for having low or undetectable levels of C-peptide (C-

peptide was detectable in 24/67 samples using the ROCHE assay and only 15/67

samples using the Mercodia assay, All patients detected by the Mercodia assay

were also detected by the Roche assay. This suggests that the Roche assay is able

to detect C-peptide at lower concentrations than the Mercodia assay.