

Supplemental Table: Primers used for quantitative PCR

L32 For	5'-AAGCGAAACTGGCGGAAAC-3'
L32 Rev	5'-TAACCGATGTTGGGCATCAG-3'
NFATc1 For	5'-GCCTCGTATCAGTGGGCGAAG-3'
NFATc1 Rev	5'-CGAAGCTCGTATGGACCA-3'
IkB α For	5'-CTGCAGGCCACCAACTACAA-3'
IkB α Rev	5'-CAGCACCCAAAGTCACCAAGT-3'
IFN γ For	5'-GTCAACAACCCACAGGTCCAG-3'
IFN γ Rev	5'-CCTTTTCCGCTTCCTGAGG-3'
IL-2 For	5'-CCTGAGCAGGATGGAGAATTACA-3'
IL-2 Rev	5'-TCCAGAACATGCCGCAGAG-3'
TNF α For	5'-AATGGCCTCCCTCTCATCAGT-3'
TNF α Rev	5'-GCTACAGGCTTGTCACCTCGAATT-3'
Cav1 α For	5'-AACATCTACAAGCCCAACAACAAGG-3'
Cav1 β For	5'-AAGCAAGTGTATGACGCGCACACC-3'
Cav1 Rev	5'-GGTTCTGCAATCACATCTTCAAAGTC-3'