

The American Journal of Human Genetics, Volume 93

Supplemental Data

Exome Sequence Reveals Mutations in CoA Synthase

as a Cause of Neurodegeneration

with Brain Iron Accumulation

Sabrina Dusi, Lorella Valletta, Tobias B. Haack, Hugo Tsuchiya, Paola Venco, Sebastiano Pasqualato, Paola Goffrini, Marco Tigano, Nikita Demchenko, Thomas Wieland, Thomas Schwarzmayr, Tim M. Strom, Federica Invernizzi, Barbara Garavaglia, Allison Gregory, Lynn Sanford, Jeffrey Hamada, Conceição Bettencourt, Henry Houlden, Luisa Chiapparini, Giovanna Zorzi, Manju A. Kurian, Nardo Nardocci, Holger Prokisch, Susan Hayflick, Ivan Gout, and Valeria Tiranti

Table S1. Next generation Sequencing Statistics

Id	Type	Reads	Mapped	Percent	Seq (Gb)	on bait	Avg cov	Cov 1x	Cov 4x	Cov 8x	Cov 20x
#55633	SureSelect50Mb	87259289	80916002	92.73	8.81	73.85	93.96	99.61	98.79	97.3	92.4

Table S2. Variants Identified by Exome Sequencing

Variants filtering	
Synonymous variants	11846
NSV	11601
NSV with frequency <0.1% in ‘in-house and public databases	340
≥ 2NSV / gene	22
Genes carrying homozygous rare NSV	12

NSV = missense, nonsense, stop/loss, splice site disruption, insertions, deletions

Table S3: detailed list of the 12 genes carrying homozygous rare NSV

gene #	gene symbol	OMIM	chromosomal position of identified variant	UCSC transcript, predicted mutation at nucleotide and protein level	dbSNP (rs number)	prefunction	pph2	Sift	Occurrence of predicted recessive-type variants in 3,159 samples with unrelated phenotypes	Other arguments against a causal role
1	GUCA2A	*139392	chr1:42629092-42629092	uc001chd.1, c.265G>T, p.Glu89*	rs147564899	nonsense	0.39	-		uc001chd.1, c.265G>T, p.Glu89* also found in healthy subjects I-2 and II-4 of family 1
2	HRNR	*114085	chr1:152189999-152189999	uc001ezt.1, c.4106G>T, p.Ser1369Phe		missense	benign	0.05	>1000 other patients with predicted recessive type mutations and unrelated phenotypes	
			chr1:152190002-152190002	uc001ezt.1, c.4103G>A, p.Gly1368Asp		missense	benign	0.3		
3	ADAM8	*602267	chr10:135082994-135082994	uc021qbe.1, c.1903T>G, p.Trp635Gly		missense	probably damaging		16 other patients with predicted recessive type mutations and unrelated phenotypes	
4	FBXO47	*609498	chr17:37118301-37118301	uc002hrc.2, c.182-1G>A, p?		splice			absence of mutations in 56 NBIA subjects	Suggested role as a tumor suppressor.
5	CACNB1	*114207	chr17:37353687-37353687	uc002hrc.2, c.62A>G, p.Glu21Gy		missense	benign	0.17	1 other patient with predicted compound heterozygous mutations and mental retardation.	
6	COASY	*609855	chr17:40717686-40717686	uc010cyj.3, c.1582C>T, p.Arg528Cys	rs140709867	missense	probably damaging	0	-	
7	BZRAP1	*610764	chr17:56385086-56385086	uc002ivx.4, c.4869G>C, p.Gln1623His	rs141269009	missense	possibly damaging	0.12	52 other patients with predicted recessive type mutations and unrelated phenotypes	
8	C17orf47	not available	chr17:56621291-56621291	uc002iwq.2, c.257C>T, p.Arg86Gln		missense	benign	0.28	2 other patients with predicted recessive type mutations and unrelated phenotypes	
9	LRP1B	*608766	chr2:141641491-141641491	uc002hv.1, c.4064A>G, p.Asn1365Ser		missense	probably damaging	0.56	151 other patients with predicted recessive type mutations and unrelated phenotypes	
10	EVC2	*607261	chr4:5624616-5624616	uc003gj.3, p.Ile171Asp		missense	possibly damaging	0.32	25 other patients with predicted recessive type mutations and unrelated phenotypes	Associated with Ellis-van Creveld syndrome (MIM *607261)
			chr4:5667370-5667370	uc003gj.3, c.877C>T, p.Pro293Ser		missense	probably damaging	0.14		
11	KIAA1797	*614606	chr9:20885216-20885216	uc003zog.1, c.2612C>T, p.Ala871Val	rs187654100	missense	benign	0.29	16 other patients with predicted recessive type mutations and unrelated phenotypes	
12	IFNW1	*147553	chr9:21141080-21141080	uc003zol.1, c.490A>G, p.Trp164Arg	rs147009804	missense	probably damaging	0	-	uc003zol.1, c.490A>G, p.Trp164Arg also found in healthy subjects I-2, II-4 and II-5 of family 1

1 **Table S4.** Primer sequences and PCR conditions to amplify the 9 exons of *COASY* gene.

Exons	Forward primer	Reverse primer	DNA polymerase	Annealing Temp
1a	GATTGGATCCCCAGCCC	CAAGAACCTCAAACGTGGCC	<i>AmpliTaq Gold®</i>	58°C
1b	TCCCATAACAGGCCTTGCC	GGCCTTAGCACCAATTGAGTG	<i>AmpliTaq Gold®</i>	60°C
2	TGCTTGCCTGTTCTTCACCT	CGGGCAAACTGAGAACCAAGTT	GoTaq® <i>Flexi</i>	59°C
3	GGAGCCTGGGTAGGAAGGG	CTCAGCCTCATGCCTGGG	GoTaq® <i>Flexi</i>	60°C
4	CCTGGCAATGCTGGAGAGTAG	GGCCTGGCTAGCTCCTCACT	GoTaq® <i>Flexi</i>	60°C
5	GGCCTCTGTGCTCAGTTGTCT	GCAAGGTGGTGTGGAAACA	GoTaq® <i>Flexi</i>	60°C
6	CCCAGAATGCCATTCCATT	GACGTTCTCCCAGGCAGAAC	GoTaq® <i>Flexi</i>	59°C
7	GAGGATGGCAACAGCTAAGGG	TTTGTCAAGAAGAATTCCAGCTC	GoTaq® <i>Flexi</i>	59°C
8	CCCCCCACACCATCCCTAC	GATGAGTCTTGGGAATGCGCT	GoTaq® <i>Flexi</i>	60°C
9	TACTGCCTGACCCTGCCCT	AAAGCTCGCCTCTGGCTAGC	<i>AmpliTaq Gold®</i>	60°C

2
3 GoTaq® *Flexi* DNA Polymerase (Promega); *AmpliTaq Gold®* (Life Technologies) were used. PCR
4 conditions were: 94°C, 1min (1 cycle); 94°C 30sec, annealing time was 30sec at the indicated
5 temperatures, 72°C 1min (30 cycles).

6

7