

SUPPLEMENTARY TABLE S1. LIST OF ANTIBODIES USED FOR IMMUNOFLUORESCENCE (IF), IMMUNOHISTOCHEMISTRY (IH), AND WESTERN BLOT (WB)

<i>Antibody</i>	<i>Host</i>	<i>Application</i>	<i>Dilution</i>	<i>Supplier</i>
B-ACTIN	Rabbit	WB	1:5,000	Santa Cruz
CD105	Mouse	IF	1:50	BD (wwwbdbiosciences.com)
CD14	Mouse	IF	1:100	Dako Cytomation (www.dako.com)
CD29	Mouse	IF	1:100	Abcam
CD34	Mouse	IF	1:100	Dako Cytomation
CD44	Mouse	IF	1:50	Abcam
CD45	Mouse	IF	1:100	Dako Cytomation
CD73	Mouse	IF	1:50	BD Biosciences
CD90	Mouse	IF	1:100	Abcam
COL1A2	Rabbit	WB	1:1,000	Abcam
HLAI	Mouse	IF	1:100	Dako Cytomation
HLAII	Mouse	IF	1:100	Dako Cytomation
OCT-4A	Mouse	IF	1:2	Santa Cruz
OCT-4A-PE	Mouse	FC	1:10	R&D Systems
SOX2	Rabbit	IF	1:100	Abcam
SOX2-PE	Mouse	FC	1:10	R&D Systems
SSEA4	Mouse	IF	1:50	Abcam
SSEA4-PE	Mouse	FC	1:10	R&D Systems
TRA-1-60	Mouse	IF, FC	1:10	Santa Cruz
VIMENTIN	Rabbit	IF, IH	1:100	Abcam
COLLAGEN X	Mouse and human	IH	1:100	Abcam
OSTEOPONTIN	Mouse and human	IH	1:100	Abcam