

Table S2 Primers used for ChIP-PCR assays

Amplicon relative to rDNA TSS (bp)	Forward primer (5' - 3')	Reverse primer (5' - 3')
+165/245	CTTGCGTGTGCTTGCTGT	GAAATCGGGAAAAACGTCTG
+6859/+6986	ACACCCGAAATACCGATACG	AGTGCGTTCGAAGTGTCGAT
+7372/+7518	CATCTGCTCTGGTCGAGGTT	GCAAGACCCAAACACACACA
+16140/+16275	CCCTTTTACACCTCCCCCTA	TTTCCAAAGGGCTACATTGC
+20690/+20857	GGAATCCGATGCACACTTTT	TGTGTGTGTGTGTGTGTGTA
+23072/+23196	AATACGCCGTGTGTGTGTGT	CAATGTGCCTTCATGATCAAA
+42801/+42909	AACTGTGCCTGTTCCCTCAC	GGCACCCAAAAACGAAAGTA
+43179/+43306	GACACAGGAGAGGGAAGTGC	CTCCCTGTACGACCTCCTTG