No. series	No. sterile plates over 30 inoculated	Cumulative QIS		Cumulative QIS		Normalized cumulative score (NCS)		Difference in NCS
		Minimum possible	Maximum possible	Inoqula	Manual	Inoqula	Manual	(Automate–Manual)
1	9	21	84	75	59	88	67	21
2	8	22	88	67	52	74	55	19
3	1	29	116	100	82	86	70	16
4	7	23	92	78	62	84	65	19
5	3	27	108	99	98	91	90	1
6	9	21	84	72	67	84	77	7
7	5	25	100	77	59	76	57	19
8	4	26	104	83	61	79	57	22
9	7	23	92	70	48	74	48	26
10	7	23	92	75	53	80	54	26
11	5	25	100	76	53	75	51	24
12	4	26	104	93	70	89	66	23
13	2	28	112	93	74	83	65	17
14	4	26	104	88	65	84	61	23
15	7	23	92	76	51	81	52	29

Table S1: Comparison of quality of inoculation between FA-Inoqula and manual method at the team level. Characteristics of the results obtained for series of 10 swabs corresponding with 30 inoculated plates. A total of 15 technicians were evaluated

NCS: Normalized cumulative score, which was calculated as follows:

(cumulative QIS (technician)-minimum QIS possible)*100/(cumulative QISmax-cumulative QISmin)