

Supplementary Table 4. Frequency of long-term culture-initiating cells (LTC-IC) in the blood of the indicated mice after G-CSF mobilization.

Treatment	Estimated LTC-IC frequency	Estimated LTC-IC per ml
Saline	1 in 26 μ l \pm 3	43.47
6OHDA	1 in 116 μ l* \pm 21	8.62

For this analysis serial dilutions of blood (30-5.9 μ l) were plated in wells (34–70 replicates per dilution) and 4-weeks later each well was assayed for the presence of CFU-C colonies. *p<0.05. LTC-IC frequency was then calculated using L-calc software (Stem Cell Technologies).