



SUPPLEMENTARY FIG. S3. Cytokine/chemokine secretions from naive CD4⁺ T cells under low or ambient O₂ culture conditions. A human proteome profiler array (Human Cytokine Array Panel A membrane; R&D Systems) was probed to assess alterations in secretions from naive CD4⁺ T cells. Data were analyzed using Graphpad Prism version 5.0 and a two-way analysis of variances was used followed by Bonferroni's post-test to evaluate statistical relevance. Data are reported as mean (μ) \pm SE of integrated density measurements using an Odyssey Li-COR scanner. * P < 0.05, ** P < 0.001, & P < 0.0001.