

SUPPLEMENTAL MATERIAL

Noman et al., <http://www.jem.org/cgi/content/full/jem.20131916/DC1>**Table S1.** Genomic oligonucleotide primers used for amplification of immunoprecipitated DNA samples from ChIP assays

Primer	Sequence
VEGF-HRE-M-FWD	CTGGGCAGCTGGCCTACCTACCTT
VEGF-HRE-M-REV	GTTAGTCAGTGGCGGGGAGTGAGA
LDHA-HRE-M-FWD	CGGGCTCCTGCTTCTGAGGCTGAG
LDHA-HRE-M-REV	CCCACGCGTGCCCAAATCTGAAC
SLC2A1-HRE-M-FWD	CCTGCCTGCGCCCTTACATC
SLC2A1-HRE-M-REV	CATACCCAAACCCGCCCTACCTA
HRE-1-PDL1-M-FWD	TCCTTATGGAACCAAGGGAAC
HRE-1-PDL1-M-REV	TGGCCTCGAACTCAGAAATC
HRE-2-3-PDL1-M-FWD	TGAGGAAGTCACCAAATCCA
HRE-2-3-PDL1-M-REV	CGATTCCCCATCTCTTTCAA
HRE-4-PDL1-M-FWD	GTTTCACAGACAGCGGAGGT
HRE-4-PDL1-M-REV	CTTTAAAGTGCCTGCAAGC