



Suppl. Fig 1

A



B



Supplementary figure 1

Neuraminidase treatment does not alter erythrocyte volume. Human erythrocytes treated with or without neuroamidase (60 mU ml^{-1}) for 60 minutes, washed, and investigated in a flow cytometer. (A) Representative experiments show side scatter (SS) as a function of electronic volume (EV) from 2500 cells during control conditions (left side image) and after neuroamidase treatment (right side image). (B) Mean values \pm SEM of side scatter and electronic volume from 4 experiments are shown in each image.