

## ADDITIONAL FILE 1

**Primers used for amplification of *CLCN5* 5' cDNA ends.** The table report the sequence of adaptor GeneRacer 5' primers supplied with the kit and the sequence of specific *CLCN5* gene antisense primers.

<i>PRIMER</i>	<i>SEQUENCE</i> (5' → 3')
GeneRacer 5' primer	CGACTGGAGCACGAGGACACTGA
r RACE Ex 2	ACACCAGGGATTGGCTCCTCCAAGAA
GeneRacer 5' nested primer	GGACACTGACATGGACTGAAGGAGTA
r RACE Ex 2 bis	CTCTCGGTGCCTATCCCGGTCTCGAG
r RACE cn 1-2	ACCATTGTACGACTTGTCCTCTGTATTG
r RACE sz 1-2	CCATTGTACGACTTGTCCTCTGGTAGG
r RACE 1c	AGGGACAAGGCTTCCTCCAGACTGTTT