

Table S2. Primers used for quantitative real time PCR.

Primer Name	Primer Sequence	Annealing Temperature (°C)
B_ACTIN_RT_F	AAAGACCTGTACGCCAACAC	59
B_ACTIN_RT_R	GTCATACTCCTGCTTGCTGAT	
CYB5R2_RT_F	GAACATGAAAATCGGGGAGA	59
CYB5R2_RT_R	TGTCTGGTTGGCAAAGATGA	
GSTP1_RT_F	ATCAGGGCCAGAGCTGGA	59
GSTP1_RT_R	ATAGGCAGGAGGCTTTGAGTGA	
SH3YL1_RT_F	GCAGGCGTGTCTTTAGAAGG	59
SH3YL1_RT_R	CTCCCTTGCTGCTTTTCTTG	
SPINT2_RT_F	AACAGCAATAATTACCTGACC	59
SPINT2_RT_R	AAGGATGCACGGCAAGGC	