

Table S2 KASP primers for marker BobWhite_c8266_227_TG_5AL. Note, the KASP design assays nucleotides on the opposing strand to the equivalent assay on the Illumina iSelect 90,000 SNP array

KASP Primer	Nucleotide assayed	Predicted phenotype	Primer sequence
Allele X	T	awnless	CGTCCATGGAGTCGTTCTCAAT
Allele Y	G	awned	GTCCATGGAGTCGTTCTCAAG
Common			GTGGTACACGTCCGGAAGAAT