

Table S2. Primers used for multiplex PCR.

| Primers | Cluster | Oligonucleotide sequence (5'-3') | Amplicon size (bp) |
|----------------|------------------|---|---------------------------|
| gyrB_fw | Cluster A | CTCTCCGCCAGGAGGTC | 466 |
| gyrB_rev | | CCGACACGGTAGGCAGGTAG | |
| narx_fw | Cluster B | TCTCGCAGAGCCTCCAGC | 76 |
| narX_rev | | CGGCTGTCGTAGGGTTCT | |
| bfmS_fw | Cluster C | ACCTGGTCCGCGAAGGCA | 100 |
| bfmS_rev | | AAGACCAGGCTGTCGAGGAA | |
| mucA_fw | Internal control | GAGTCGTGAAGCCCTGCA | 333 |
| mucA_rev | | AGGGCGTCGTTCTGGTTGTA | |