



Figure S1. Amplicon with or without homologous M13 adaptor sequence strategy. The amplicon sequence alignments are shown demonstrating there is no homology among them; nucleotides are colored red, blue, yellow, and green for A, C, G, and T, respectively. PCR products with or without M13 adaptors for potato blight RXLR, giant squid ND4, pigeon COI and grape rbcL were captured separately with either biotinylated (red circle) 100 bp PCR products for pigeon COI or blight RXLR. Resulting capture libraries were quantified by qPCR for each amplicon individually.