

Table S3. Pyrosequencing primers for *ASIP* analysis.

<i>ASIP</i> Variant	Amplicon	
	Size (bp)	Primers
Indel and A127G	105	F: Biotin TGAGGAAAAACCCAGAGATGAC R: GGGGTCAAGCTGGGCTACTTA Seq: ATAGAGACAGAAGRGAGATCCAACAT
T142C	105	F: Biotin TGAGGAAAAACCCAGAGATGAC R: GGGGTCAAGCTGGGCTACTTA Seq: GGGCTACTTACCTACWATAGAGACAGAAG
G41C	62	F: Biotin TGGCCACCCTGCTGGTCT R: TTCCTCAGGTGCCAGGTGAC Seq: GTAGGCAGTGAGGAGGCACAG