

Table S5. Primers and adapters used for generation of AFLP markers.

Primers/adapters	Sequence (5'--> 3')	5'-Mod.
<i>Adapters</i>		
Apal_fwd	TCGTAGACTGCGTACAGGCC	
Apal_rev	TGTACGCAGTCTAC	
TaqI_fwd	GACGATGAGTCCTGAC	
TaqI_rev	CGGTCAGGACTCAT	
<i>Preamplification primers</i>		
Apal	GACTGCGTACAGGCC	
TaqI	GATGAGTCCTGACCGA	
<i>Selective amplification primers</i>		
TaqI-G_IR700	GATGAGTCCTGACCGAG	IRD700
TaqI-C_IR800	GATGAGTCCTGACCGAC	IRD800
Apal-CA	GACTGCGTACAGGCCCA	
Apal-CC	GACTGCGTACAGGCCCC	
Apal-CG	GACTGCGTACAGGCCCCG	
Apal-CT	GACTGCGTACAGGCCCT	
Apal-GA	GACTGCGTACAGGCCGA	
Apal-GC	GACTGCGTACAGGCCGC	
Apal-GG	GACTGCGTACAGGCCGG	
Apal-GT	GACTGCGTACAGGCCGT	
Apal-AA	GACTGCGTACAGGCCAA	
Apal-AC	GACTGCGTACAGGCCAC	