

Supplemental Table 1: Primers used for qRT-PCR

Amplicon	Forward Primer (5'-3')	Reverse Primer (5'-3')
<i>mHoxa7</i>	CACATGGTTCCAGTTTGTGG	CTGAACCCCAGATCCACTGT
<i>mHoxa9</i>	CCCCGACTTCAGTCCTTGC	GATGCACGTAGGGGTGGTG
<i>mHoxa10</i>	GGAAGCATGGACATTCAGGT	CCAGGCAAGCAAGACCTTAG
<i>mHoxa11</i>	CCCCTGGTGGTTCACTCTTA	CTTGGGGCACACAGTTTCTT
<i>mGapdh</i>	CCTGGAGAAACCTGCCAAGTATG	AGAGTGGGAGTTGCTGTTGAAGTC
<i>mβ2M</i>	ACCGGCCTGTATGCTATCCAGAAA	GGTGAATTCAGTGTGAGCCAGGAT
<i>hHOXA9</i>	CTGTCCCACGCTTGACACTC	CTCCGCCGCTCTCATTCTC
<i>hHOXA10</i>	GCCCCTCCGAGAGCAGCAAAG	AGGTGGACGCTGCGGCTAATCTCTA
<i>hGAPDH</i>	ATGGGGAAGGTGAAGGTCG	GGGGTCATTGATGGCAACAATA
<i>hβ2M</i>	GTGTCTGGGTTTCATCCATC	GGCAGGCATACTCATCTTTT