

Table S2. Primers used in Quantitative PCR

FLIP _L forward	5'-TGCCTCAGAGCATACCTGAAGAGA-3'
FLIP _L reverse	5'-TCTCGAAGAAGCTCTGTCTCATTGC-3'
GAPDH forward	5'-ATGACATCAAGAAGGTGGTGA-3'
GAPDH reverse	5'-CTGTAGCCAAATTCGTTGTCA-3'
IL6 receptor forward	5'-GGGCTGGAACGGTCAAAGACA-3'
IL6 receptor reverse	5'-GGATGACACAGTGATGCTGGA-3'
36B4 forward	5'-GATTGGCTACCCAAGTGTG-3'
36B4 reverse	5'-CAGGGGCAGCAGCCACAAA-3'
IL6 forward	5'-CCAGGAGCCCAGCTATGAAC-3'
IL6 reverse	5'-CCCAGGGAGAAGGCAACTG-3'
IL8 forward	5'-TTGGCAGCCTTCCTGATTTC-3'
IL8 reverse	5'-TCTTTAGCACTCCTTGGCAAAC-3'