


(A)

Primers used for PCR amplification of the *rps4* fragment insert

Primer name	Sequence	Location	Primer combination	Expected product size
Forward 1	5' -TGCTGGAATTGTATTCCT-3'	29,198-29,215	F1-R1	510 bp
Forward 2	5' -TCCTCCTACTAATCTAATAG-3'	29,431-29,450	F2-R1	277 bp
Reverse 1	5' -ACGGTATTTTAACTTAC-3'	29,691-29,708	F1-R3	635 bp
Reverse 3	5' -GAGCCACTGATTTGTAATCAG-3'	29,813-29,833	F2-R3	402 bp

(B)



(C)

