

Table S1. Intra-day assaying of lyso-Gb3 in charcoal-treated plasma.

Spiked concentration (nM)	Determined concentration (nM)	Mean determined concentration (nM)	S.D. (nM)	Precision (%)	Accuracy (%)
0.080	0.090 0.082 0.088 0.086 0.090	0.087	0.003	3.7	9.1
0.40	0.35 0.36 0.37 0.38 0.36	0.36	0.01	3.0	-9.7
10	9.0 8.7 9.3 9.1 9.2	9.1	0.2	2.4	-9.5
200	205 198 206 206 204	204	3	1.6	2.0