

**S1 Table: Primary and secondary antibodies used for the immunohistological stainings**

Primary Antibody	Source	Dilution
Polyclonal rabbit $\alpha$ -SMA	Abcam	1:100
Polyclonal goat anti-collagen I	Southern Biotech	1:100
Polyclonal goat anti-collagen III	Southern Biotech	1:100
Polyclonal rabbit Cytokeratin 19	Novus Biologicals	1:100
Polyclonal goat anti-desmin	Santa Cruz	1:100
Polyclonal rat anti-F4/80 (MCA497)	AbD Serotec	1:500
Polyclonal rabbit anti-Ki67	Millipore	1:100
Secondary Antibody	Source	Dilution
Polyclonal Goat anti-rabbit IgG	DAKO	1:100
Polyclonal Rabbit anti-rat IgG	DAKO	1:100
Polyclonal Rabbit anti-Goat IgG	DAKO	1:100