



Supplementary information, Figure S11 Test for fluorescence bleed-through between different channels by single-labeled specimens. HeLa cells were transfected with caveolin-1-CFP (**A**), caveolin-1-GFP (**B**), or caveolin-1-mRFP (**C**) and imaged by the Olympus FluoView FV1000-IX81 confocal microscope. For the four color imaging, the lasers and the corresponding emission filters for each channel used were shown at the top of the images. The fluorescence bleed-through between different channels is negligible under our experiment condition. Scale bar, 20 μm .