

Supplementary Table 1. Primers used for real-time PCR analyses

<b>Gene</b>		<b>Primer</b>
mCcmd1	Forward	5'-GCGTACCCTGACACCAATCT-3'
	Reverse	5'-ATCTCCTTCTGCACGCACTT-3'
mCcmd2	Forward	5'-TCGATGATTGCAACTGGAAG-3'
	Reverse	5'-AGAGCTTCGATTTGCTCCTG-3'
mE2f3	Forward	5'-TCTAAAGACCCCAAGGGCA-3'
	Reverse	5'-GATCCAGGACCCCATCAGGA-3'
mRbl1	Forward	5'-CCCGGAGCCGAAAACAGAG-3'
	Reverse	5'-ACCAATCATACGGAAATCCCC-3'
mCcnb2	Forward	5'-TGTCACAAGCAGCCGAAAC-3'
	Reverse	5'-TCAGAGAAAGCTTGGCAGAGG-3'
mCdk1	Forward	5'-AGAAGGTACTTACGGTGTGGT-3'
	Reverse	5'-GAGAGATTTCCCGAATTGCAGT-3'
mChek1	Forward	5'-GTTAAGCCACGAGAATGTAGTGA-3'
	Reverse	5'-TTCAGGCATCCCTATGTCTGG-3'
mCks1b	Forward	5'-TATTTCGGACAAATACGACGACG-3'
	Reverse	5'-GGTTCCTCCATTTCAGATTCAGAC-3'
mCks2	Forward	5'-TCGATGAGCACTACGAGTACC-3'
	Reverse	5'-CCATCCTAGACTCTGTTGGACAC-3'
mRprm	Forward	5'-GTGTGGTGCAGATCGCAGT-3'
	Reverse	5'-ATCATGCCTTCGGACTTGATG-3'
mTop2a	Forward	5'-CAACTGGAACATATACTGCTCCG-3'
	Reverse	5'-GGGTCCCTTTGTTTGTATCAGC-3'
mCcna2	Forward	5'-GCAGTTTTGAATCACCACATGC-3'
	Reverse	5'-TGGCTGCCTCTTCATGTAACC-3'
mBmp2	Forward	5'-AGATCTGTACCGCAGGCACT-3'
	Reverse	5'-CGTCACTGGGGACAGAACTT-3'
m18s	Forward	5'-GTAACCCGTTGAACCCATT-3'
	Reverse	5'-CCATCCAATCGGTAGTAGCG-3'
mFkbp5		Mm00487401
mI13ra2		Mm00515166
mCited2		Mm00516121
m18s		4319413E