

Primers used:

Forward primer:

All bacteria Amp_27F_v4 5'-AGRGTTYGATYMTGGCTCAG-3'

Reverse primers:

All bacteria Eub_338R_3W 5'-GCWGCCWCCCGTAGGWGT-3'

Bifidobacterium Bif_162R 5'-CCGGYATTACCACCCGTTT-3'

Collinsella Coll_230R 5'-TCACCCCGCCGTCTACCTGAT-3'

Faecalibacterium Faec_396R 5'-CCGAAGACCTTCTTCCTCC-3'