

SUPPLEMENTAL TABLE 1. FilmArray® BioThreat Panel and qRT-PCR testing results for whole blood, plasma and urine from patients with EVD.<sup>a,b</sup>

Specimen		Result		
FilmArray®	qRT-PCR	FilmArray®	NP2 qRT-PCR	CT
Whole blood	Whole blood	Detected	Positive	20
		Detected	Positive	21
		Detected	Positive	22
		Detected	Positive	23
		Detected	Positive	24
		Detected	Positive	35
		Detected	Positive	35
		Detected	Positive	36
		Detected	Positive	36
		Detected	Positive	36
		Detected	Positive	36
		Detected	Positive	36
		Detected	Positive	37
		Detected	Positive	37
		Detected	Positive	37
		Detected	Positive	37
		Detected	Positive	37
		Not Detected	Positive	38
		Detected	Positive	39
		Not Detected	Positive	39
		Not Detected	Positive	39
		Detected	Positive	40
		Detected	Negative	NA
Whole blood	Plasma	Not Detected	Negative	NA
		Not Detected	Negative	NA
		Not Detected	Negative	NA
		Not Detected	Negative	NA
		Not Detected	Negative	NA
		Detected	Positive	23
		Detected	Positive	26
		Detected	Positive	27
		Detected	Positive	28
		Detected	Positive	31
		Detected	Positive	32
		Detected	Positive	36
		Detected	Positive	37
		Detected	Positive	38
		Detected	Positive	38
		Detected	Positive	39
		Detected	Positive	39
		Detected	Positive	39
		Detected	Negative	NA
		Detected	Negative	NA
Urine	Urine	Not Detected	Negative	NA
		Not Detected	Negative	NA
		Not Detected	Negative	NA
		Not Detected	Negative	NA
		Not Detected	Negative	NA

Detected	Positive	30
Detected	Positive	34
Detected	Positive	35
Detected	Positive	37
Detected	Positive	38
Detected	Positive	38
<b>Not Detected</b>	<b>Positive</b>	<b>39</b>
Detected	Positive	39
<b>Detected</b>	<b>Negative</b>	NA
Not Detected	Negative	NA
Not Detected	Negative	NA
Not Detected	Negative	NA

Abbreviation: NA, not applicable

<sup>a</sup>Testing was performed using both the FilmArray® BioThreat panel and the CDC real-time reverse transcriptase polymerase chain reaction (qRT-PCR) assay.

<sup>b</sup>Discrepancies in bold.