

Preservation of prostate tissue in fidelity

Supplementary Table 1. Sequences of primers used for real time PCR assay

| Name | Sequence (5'-3') |
|------------------|------------------------|
| GAPDH-F | AGTCCCAGAAACAGGAGGT |
| GAPDH-R-88 bp | AGAGCGCGAAAGGAAAGAA |
| GAPDH-R-270 bp | GTTAGTCACCGGCAGGCTTT |
| GAPDH-R-420 bp | CACAAGAGGACCTCCATAAACC |
| GAPDH-R-486 bp | ACCCATGACTCAGCTTCTCC |
| GAPDH-R-693 bp | GTTTCCGGAAGACGGAATG |
| GAPDH-R-817 bp | GGGAGCACAGGTAAGTGCAT |
| B-ACTIN-F | AAATCTGGCACCACACCTTC |
| B-ACTIN-R-95 bp | GTTGGCCTTGGGGTTTCAG |
| B-ACTIN-R-626 bp | AGAGGCGTACAGGGATAGCA |
| B-ACTIN-R-703 bp | AGGGCATACCCTCGTAGAT |
| B-ACTIN-R-828 bp | TAATGTCACGCACGATTTC |
| B-ACTIN-R-983 bp | AGGCTGGAAGAGTGCCTCAG |