

S2 Table. Quantitative-PCR primers

Amplicon	Forward primer (5'-3')	Reverse primer (5'-3')	Probe (5'-3')
Strong stop	TAACTAGGGAACCCACTGC	GCTAGAGATTTTCCACACTG	FAM- ACACAACAGACGGGCACAC ACTA-TAMRA
Second strand	CTGCGTCGAGAGAGCTCCTC TGGTT	TAACTAGGGAACCCACTGC	FAM- ACACAACAGACGGGCACAC ACTA-TAMRA