

ENZO Testosterone Plate Layout

A1 Blank	A2 Std 5	A3	A4	A5	A6	A7	A8	A9	A10	A11	A12 Std 4
B1 TA	B2 Std 6	B3	B4	B5	B6	B7	B8	B9	B10	B11	B12 Std 3
C1 NSB	C2 Std 7	C3	C4	C5	C6	C7	C8	C9	C10	C11 Low Control	C12 Std 2
D1 Bo	D2 Low Control	D3	D4	D5	D6	D7	D8	D9	D10	D11 Medium Control	D12 Std 1
E1 Std 1	E2 Medium Control	E3	E4	E5	E6	E7	E8	E9	E10	E11 High Control	E12 Bo
F1 Std 2	F2 High Control	F3	F4	F5	F6	F7	F8	F9	F10	F11 Std 7	F12 NSB
G1 Std 3	G2	G3	G4	G5	G6	G7	G8	G9	G10	G11 Std 6	G12 TA
H1 Std 4	H2	H3	H4	H5	H6	H7	H8	H9	H10	H11 Std 5	H12 Blank

Kit Lot #: _____ Exp. Date: _____ Date Run: _____ Tech: _____

1st Incubation: Start time: _____

Notes: _____

2nd Incubation: Start time: _____

End time: _____

3rd Incubation: Start time: _____

End time: _____

Pipette Calib 100 uL mili-Q H₂O = _____ g Pipette: _____

Pipette Calib 40 uL mili-Q H₂O = _____ g Pipette: _____