

S10 Table. LNAs used for miRNA in situ hybridization and primers used for probe preparation

		size
miR-31	/5DigN/AGCTTCGCCGTAATCTTGCCA/3Dig_N/	
Bantam	/5DigN/ACCAGCTTTAATCGCGATCTCA/3Dig_N/	
Scrambled Sj miRNA	/5DigN/GTGTAACACGTCTATACGCCCA/3Dig_N/	
Frizz7 (EU370927)	Forward primer: 5'TGTTGGCGAAGCATTCAACG3` Reverse primer: 5` ATCGAAATTAATACGACTCACTATAGGGCCGACTGAGAAGGTGAGCGAA3`	661nt
Smad1 (AY815078.1)	Forward primer: 5`TGGGGCATGCAACTCAAGAT3` Reverse primer: 5` ATCGAAATTAATACGACTCACTATAGGGC CTGAGAGACTGTGGTGCTGG3`	336nt