

## Mismatch repair gene defects in sporadic colorectal cancer enhance immune surveillance

### Supplementary Material

**Supplementary Table 1.** Features of the antibodies used for IHC

	firm	clone	dilution	unmasking
CD80	RD System		1:100	pH 6
CD4	Dako	4B12	Ready to use	pH 9
CD8	Thermo Scientific	sp16	1:50	pH 9
T-bet	Santa Cruz	H-210	1:50	pH 9
Fox-P3	Abcam		1:100	pH 9
MLH1	Dako	ES05	READY TO USE	pH 9
MSH2	Thermo Scientific	25D12	READY TO USE	pH 9
MSH2	Dako	FE11	READY TO USE	pH 9
MSH6	Dako	EP49	1:50	pH 9
MSH6	Dako	EP51	1:50	pH 9