

S2 Table. Probe sequences for TaqMan SNP assays

Variant	Forward primer (5'-3')	Reverse primer (3'-5')	Reporter 1 (5'-3')	Reporter 2 (5'-3')
c.3010-1G>C	CCAGGCAATCCAATAAGTAAAAACATACAT	GGCTCCCCATCATCAAAAGCA	VIC-CCACCAACTTAAAAAG-NFQ	FAM-CCACCAAGTAAAAAG-NFQ
p.D642N	TCAATTACTTGATAATTGAATTATGTTTTCTCTC	CAATGGCATATGTTATTGAGTCTCCCT	VIC-TGAAGGCAACTGATCGAG-NFQ	FAM-AAGGCAACTAATCGAG-NFQ
p.V469A	GACCCAGAGCTTCACCTTTTTCT	GAAAAGGTGTAAGTTTGCTGTTCTT	VIC-CCCTGTCCACTGGTTG-NFQ	FAM-CCTGTCCGCTGGTTG-NFQ
p.T281A	TGTGCCAAACACTCGTGATTG	CAACTCAGGTATGGCAGCTTG	VIC-CCGTCCACTCACTTAT-NFQ	FAM-CGTCCACTCGCTTAT-NFQ
p.R219K	CTTGGATTATGACAAAGTAGCGAGTCT	ACCTTTGAAATCCCCTAATGTTGACT	VIC-CATAGTTGAGCCTCTTCC-NFQ	FAM-ATAGTTGAGCTTCTTCC-NFQ
p.M60I	AGTAGAATATGTTTTCAAGTACCTTTCTTTGG	CAGGATCCATCAACACCCAGTAATC	VIC-TGGACAACATGCTGATCA-NFQ	FAM-TGGACAACATACTGATCA-NFQ