

Table S2. Primers used for quantitative qRT-PCR.

Gene	Forward primer (5' > 3')	Reverse primer (5' > 3')
<i>Trib1</i>	CTTCAAGCAGATTGTCTCCGC	CTAAGCTGGGTTCTCTCCTCC
<i>Trib2</i>	AGCCAGACTGTTCTACCAGA	GGCGTCTTCCAGGCTTTCCA
<i>Abl</i>	TGGAGATAACACTCTAAGCATAACTAAAGGT	GATGTAGTTGCTTGGGACCCA
<i>B2m</i>	TTGTCTTTCAGCAAGGACTGG	ATGCGGCATCTTCTAACCTCC
<i>Enox2</i>	GAGCTGGAGGGAACCTGATTT	CACTGGCACTACCAAACCTGCA
<i>Rnf20</i>	GGTGTCTCTTCAACGGAGGAA	TAGTGAGGCATCATCAGTGGC