

Table S3. Primers used for detection of endogenous *Tcrb* rearrangements.

	<i>Tcrb</i> locus	5'>3' sequence	Annealing temperature (°C)
Forward primer	Dβ2	GTAGGCACCTGTGGGGAAGAACT	58
	Vβ2	GGGTCACTGATACGGAGCTG	58
	Vβ4	GGACAATCAGACTGCCTCAAGT	58
	Vβ5.1	GTCCAACAGTTTGATGACTATCAC	56
	Vβ8	GATGACATCATCAGGTTTTGTC	56
	Vβ14	CTTCTACCTCTGTGCCTGGAGT	58
Reverse primer	Jβ2	TGAGAGCTGTCTCCTACTATCGATT	according to forward primer